

For Research Use Only

ZBTB7A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83616-4-PBS

Basic Information

Catalog Number:

83616-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC113511

GeneID (NCBI):

51341

UNIPROT ID:

O95365

Full Name:

zinc finger and BTB domain containing 7A

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

240555H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00594-1, ZBTB7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83616-1-PBS. Detection antibody: 83616-4-PBS. Standard: SY00194. Range: 39-5000 pg/mL



Cytometric bead array standard curve of MP00594-2, ZBTB7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83616-2-PBS. Detection antibody: 83616-4-PBS. Standard: Ag38535. Range: 0.313-40 ng/mL