

For Research Use Only

# SLC35F3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83617-1-PBS

## Basic Information

Catalog Number: 83617-1-PBS	GenBank Accession Number: NM_001300845	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 148641	CloneNo.: 240386E2
Source: Rabbit	UNIPROT ID: Q8IY50	
Isotype: IgG	Full Name: solute carrier family 35, member F3	
Immunogen Catalog Number: AG35308	Calculated MW: 421 aa, 47 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00596-1, SLC35F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83617-1-PBS. Detection antibody: 83617-2-PBS. Standard: Ag35308. Range: 0.78-100 ng/mL.