

For Research Use Only

# RBMS1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83623-1-PBS



## Basic Information

Catalog Number:

83623-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1493

GenBank Accession Number:

BC018951

GeneID (NCBI):

5937

UNIPROT ID:

P29558

Full Name:

RNA binding motif, single stranded  
interacting protein 1

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

240677E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00575-1, RBMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83623-4-PBS. Detection antibody: 83623-1-PBS. Standard: Ag1493. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00575-2, RBMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83623-3-PBS. Detection antibody: 83623-1-PBS. Standard: Ag1493. Range: 0.78-100 ng/mL