

For Research Use Only

RBMS1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83623-3-PBS



Basic Information

Catalog Number:

83623-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1493

GenBank Accession Number:

BC018951

GeneID (NCBI):

5937

UNIPROT ID:

P29558

Full Name:

RNA binding motif, single stranded
interacting protein 1

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

240677F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00575-2, RBMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83623-3-PBS. Detection antibody: 83623-1-PBS. Standard: Ag1493. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00575-3, RBMS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83623-4-PBS. Detection antibody: 83623-1-PBS. Standard: Ag1493. Range: 0.78-100 ng/mL