

For Research Use Only

PD-1/CD279 Recombinant antibody, PBS Only (Capture)

Catalog Number:83632-2-PBS



Basic Information

Catalog Number: 83632-2-PBS	GenBank Accession Number: BC074740	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5133	CloneNo.: 240724G7
Source: Rabbit	UNIPROT ID: Q15116	
Isotype: IgG	Full Name: programmed cell death 1	
	Calculated MW: 288 aa, 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00612-1, PD-1/CD279 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83632-2-PBS. Detection antibody: 83632-1-PBS. Standard: Eg0974. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00612-2, PD-1/CD279 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83632-2-PBS. Detection antibody: 83632-4-PBS. Standard: Eg0974. Range: 0.313-40 ng/mL.