

For Research Use Only

ATP6V1E1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83636-5-PBS

Basic Information

Catalog Number:

83636-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7399

GenBank Accession Number:

BC004443

GeneID (NCBI):

529

UNIPROT ID:

P36543

Full Name:

ATPase, H⁺ transporting, lysosomal
31kDa, V1 subunit E1

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

240584F7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00628-3, ATP6V1E1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83636-5-PBS. Detection antibody: 83636-3-PBS. Standard: Ag7399. Range: 2.5-40 ng/mL