

For Research Use Only

RCN1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83638-2-PBS



Basic Information

Catalog Number:	83638-2-PBS	GenBank Accession Number:	BC010120	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	5954	CloneNo.:	240473H1
Source:	Rabbit	UNIPROT ID:	Q15293		
Isotype:	IgG	Full Name:	reticulocalbin 1, EF-hand calcium binding domain		
Immunogen Catalog Number:	AG13256	Calculated MW:	39 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00637-1, RCN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83638-2-PBS. Detection antibody: 83638-1-PBS. Standard: Ag13256. Range: 0.313-40 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83638-2-PBS against Human RCN1 was performed. The affinity constant is 1.26 nM.