

For Research Use Only

# WFS1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83642-2-PBS

## Basic Information

|                                      |  |  |
|--------------------------------------|--|--|
| Catalog Number:<br>83642-2-PBS       | GenBank Accession Number:<br>BC030130        | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>7466                       | CloneNo.:<br>240538B10                         |
| Source:<br>Rabbit                    | UNIPROT ID:<br>O76024                        |  |
| Isotype:<br>IgG                      | Full Name:<br>Wolfram syndrome 1 (wolframin) |  |
| Immunogen Catalog Number:<br>AG25735 | Calculated MW:<br>890 aa, 100 kDa            |  |

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00629-2, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-2-PBS. Detection antibody: 83642-4-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00629-3, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-2-PBS. Detection antibody: 83642-3-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL