

For Research Use Only

# WFS1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83642-3-PBS

## Basic Information

Catalog Number: 83642-3-PBS	GenBank Accession Number: BC030130	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7466	CloneNo.: 240538B2
Source: Rabbit	UNIPROT ID: O76024	
Isotype: IgG	Full Name: Wolfram syndrome 1 (wolframin)	
Immunogen Catalog Number: AG25735	Calculated MW: 890 aa, 100 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00629-1, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-1-PBS. Detection antibody: 83642-3-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00629-3, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-2-PBS. Detection antibody: 83642-3-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL