

For Research Use Only

# HEATR3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83643-3-PBS

## Basic Information

Catalog Number: 83643-3-PBS	GenBank Accession Number: BC033077	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55027	CloneNo.: 240561A12
Source: Rabbit	UNIPROT ID: Q7Z4Q2	
Isotype: IgG	Full Name: HEAT repeat containing 3	
Immunogen Catalog Number: AG26328	Calculated MW: 680 aa, 75 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00631-2, HEATR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83643-1-PBS. Detection antibody: 83643-3-PBS. Standard: Ag26328. Range: 0.391-50 ng/mL.