

For Research Use Only

HEATR3 Recombinant antibody, PBS Only

Catalog Number: 83643-4-PBS

Featured Product



Basic Information

Catalog Number:

83643-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26328

GenBank Accession Number:

BC033077

GeneID (NCBI):

55027

UNIPROT ID:

Q7Z4Q2

Full Name:

HEAT repeat containing 3

Calculated MW:

680 aa, 75 kDa

Observed MW:

70-75 kDa

Purification Method:

Protein A purification

CloneNo.:

240561F7

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83643-4-RR (HEATR3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83643-4-PBS in a different storage buffer formulation.



WB result of HEATR3 antibody (83643-4-RR; 1:400; incubated at room temperature for 1.5 hours) with sh-Control and sh-HEATR3 transfected HEK-293 cells. This data was developed using the same antibody clone with 83643-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83643-4-RR against Human HEATR3 were performed. The affinity constant is 0.657 nM.