

For Research Use Only

Mouse IL-1 alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 83644-3-PBS



Basic Information

Catalog Number:

83644-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_010554

GeneID (NCBI):

16175

UNIPROT ID:

P01582

Full Name:

interleukin 1 alpha

Calculated MW:

31kd

Purification Method:

Protein A purification

CloneNo.:

240431B10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00634-1, MOUSE IL-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83644-4-PBS. Detection antibody: 83644-3-PBS. Standard: Eg0263. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00634-3, MOUSE IL-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83644-2-PBS. Detection antibody: 83644-3-PBS. Standard: Eg0263. Range: 0.78-50 ng/mL