

For Research Use Only

# Mouse IL-1 alpha Recombinant antibody, PBS Only (Capture)

Catalog Number: 83644-4-PBS



## Basic Information

Catalog Number:

83644-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_010554

GeneID (NCBI):

16175

UNIPROT ID:

P01582

Full Name:

interleukin 1 alpha

Calculated MW:

31kd

Purification Method:

Protein A purification

CloneNo.:

240431G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00634-1, MOUSE IL-1 alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83644-4-PBS. Detection antibody: 83644-3-PBS. Standard: Eg0263. Range: 0.78-100 ng/mL