

For Research Use Only

PLA2G4E Recombinant antibody, PBS Only (Capture)



Catalog Number: 83645-1-PBS

Basic Information

Catalog Number:	83645-1-PBS	GenBank Accession Number:	BC101584	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	123745	CloneNo.:	240620B4
Source:	Rabbit	UNIPROT ID:	Q3MJ16		
Isotype:	IgG	Full Name:	phospholipase A2, group IVE		
Immunogen Catalog Number:	AG12808	Calculated MW:	868 aa, 99 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00638-2, PLA2G4E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83645-1-PBS. Detection antibody: 83645-3-PBS. Standard: Ag12808. Range: 0.625-80 ng/mL.