

For Research Use Only

# SLC5A8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83652-1-PBS

## Basic Information

Catalog Number: 83652-1-PBS	GenBank Accession Number: BC110492	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 160728	CloneNo.: 240592B3
Source: Rabbit	UNIPROT ID: Q8N695	
Isotype: IgG	Full Name: solute carrier family 5 (iodide transporter), member 8	
Immunogen Catalog Number: AG30183	Calculated MW: 610 aa, 67 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00625-1, SLC5A8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83652-1-PBS. Detection antibody: 83652-2-PBS. Standard: Ag30183. Range: 0.78-100 ng/mL