

For Research Use Only

# CEA/CD66e Recombinant antibody, PBS Only

Catalog Number: 83653-2-PBS



## Basic Information

Catalog Number:

83653-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0449

GenBank Accession Number:

NM\_004363

GeneID (NCBI):

1048

UNIPROT ID:

P06731-1

Full Name:

carcinoembryonic antigen-related  
cell adhesion molecule 5

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

240666E11

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83653-2-PBS against Human CEA/CD66e were performed. The affinity constant is 0.910 nM.