

For Research Use Only

CEA/CD66e Recombinant antibody, PBS Only (Capture)

Catalog Number: 83653-4-PBS



Basic Information

Catalog Number:

83653-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004363

GeneID (NCBI):

1048

UNIPROT ID:

P06731-1

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 5

Calculated MW:

77 kDa

Purification Method:

Protein A purification

CloneNo.:

240666F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00626-1, CEA/CD66e Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83653-4-PBS. Detection antibody: 83653-3-PBS. Standard: Eg0449. Range: 0.625-80 ng/mL.