

For Research Use Only

LPXN Recombinant antibody, PBS Only



Catalog Number: 83657-5-PBS

Basic Information

Catalog Number:

83657-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34956

GenBank Accession Number:

BC019035

GeneID (NCBI):

9404

UNIPROT ID:

O60711

Full Name:

leupaxin

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

240483G10

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83657-5-RR (LPXN antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83657-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83657-5-RR against Human LPXN were performed. The affinity constant is 0.187 nM.