

For Research Use Only

Glucagon Recombinant antibody, PBS Only (Detector)

Catalog Number: 83659-4-PBS



Basic Information

Catalog Number:

83659-4-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8677

GenBank Accession Number:

BC005278

GeneID (NCBI):

2641

UNIPROT ID:

P01275

Full Name:

glucagon

Calculated MW:

180 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

240595C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00614-2, Glucagon Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83659-5-PBS. Detection antibody: 83659-4-PBS. Standard: Ag8677. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00614-3, Glucagon Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83659-1-PBS. Detection antibody: 83659-4-PBS. Standard: Ag8677. Range: 0.313-40 ng/mL