

For Research Use Only

CDK9 Recombinant antibody, PBS Only

Catalog Number: 83662-1-PBS



Basic Information

Catalog Number:

83662-1-PBS

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2318

GenBank Accession Number:

BC001968

GeneID (NCBI):

1025

UNIPROT ID:

P50750

Full Name:

cyclin-dependent kinase 9

Calculated MW:

372 aa, 43 kDa

Purification Method:

Protein A purification

CloneNo.:

240600A10

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assay of 83662-1-PBS against Human CDK9 was performed. The affinity constant is below 1 pM.