

For Research Use Only

FTCD Recombinant antibody, PBS Only (Detector)

Catalog Number: 83664-3-PBS



Basic Information

Catalog Number:

83664-3-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16646

GenBank Accession Number:

BC136395

GeneID (NCBI):

10841

UNIPROT ID:

O95954

Full Name:

formiminotransferase

cyclodeaminase

Calculated MW:

541 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

240654A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00636-1, FTCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83664-2-PBS. Detection antibody: 83664-3-PBS. Standard: Ag16646. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00636-2, FTCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83664-1-PBS. Detection antibody: 83664-3-PBS. Standard: Ag16646. Range: 5-80 ng/mL