

For Research Use Only

FTCD Recombinant antibody, PBS Only (Detector)

Catalog Number: 83664-4-PBS



Basic Information

Catalog Number:

83664-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16646

GenBank Accession Number:

BC136395

GeneID (NCBI):

10841

UNIPROT ID:

O95954

Full Name:

formiminotransferase

cyclodeaminase

Calculated MW:

541 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

240654C9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00636-3, Human FTCD Recombinant Matched Antibody Pair - PBS only. 83664-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag16646. 83664-4-PBS was HRP conjugated as the detection antibody. Range: 125-8000 pg/mL



The mean FTCD concentration was determined to be 114.0 pg/mL in HepG2 cell extract based on a 3.0 mg/mL extract load.