

For Research Use Only

SLITRK6 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 83665-6-PBS



Basic Information

Catalog Number:

83665-6-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5865

GenBank Accession Number:

NM_032229.2

GeneID (NCBI):

84189

UNIPROT ID:

Q9H5Y7

Full Name:

SLIT and NTRK-like family, member 6

Calculated MW:

95 kDa

Purification Method:

Protein A purification

CloneNo.:

261536C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SLITRK6 concentration was determined to be 3.7 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, THP-1 cell extract were served as a negative control.



Cytometric bead array standard curve of MP00640-2, SLITRK6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83665-6-PBS. Detection antibody: 83665-5-PBS. Standard: Eg5865. Range: 0.156-20 ng/mL.