

For Research Use Only

SLC4A7 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83676-4-PBS

Basic Information

Catalog Number:

83676-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30545

GenBank Accession Number:

NM_003615

GeneID (NCBI):

9497

UNIPROT ID:

Q9Y6M7

Full Name:

solute carrier family 4, sodium bicarbonate cotransporter, member 7

Calculated MW:

136 kDa

Purification Method:

Protein A purification

CloneNo.:

240689G11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00654-2, SLC4A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83676-1-PBS. Detection antibody: 83676-4-PBS. Standard: Ag30545. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00654-3, SLC4A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83676-3-PBS. Detection antibody: 83676-4-PBS. Standard: Ag30545. Range: 0.313-40 ng/mL