

For Research Use Only

# PBRM1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83677-5-PBS



## Basic Information

Catalog Number:

83677-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33623

GenBank Accession Number:

BC015323

GeneID (NCBI):

55193

UNIPROT ID:

Q86U86

Full Name:

polybromo 1

Calculated MW:

180 kDa

Purification Method:

Protein A purification

CloneNo.:

240655F12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00656-2, PBRM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83677-5-PBS. Detection antibody: 83677-1-PBS. Standard: Ag33623. Range: 0.625-80 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83677-5-PBS against Human PBRM1 was performed. The affinity constant is below 1 pM.