

For Research Use Only

# FSTL3/FLRG Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83681-4-PBS



## Basic Information

Catalog Number:

83681-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0949

GenBank Accession Number:

NM\_005860.3

GeneID (NCBI):

10272

UNIPROT ID:

O95633-1

Full Name:

folistatin-like 3 (secreted glycoprotein)

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

261558A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean FSTL3/FLRG concentration was determined to be 15.8 ng/mL in human serum 1, 4.4 ng/mL in human serum 2, 11.3 ng/mL in human serum 3.



Cytometric bead array standard curve of MP00667-2, FSTL3/FLRG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83681-5-PBS. Detection antibody: 83681-4-PBS. Standard: Eg0949. Range: 0.156-10 ng/mL.