

For Research Use Only

CDCA7L Recombinant antibody, PBS Only (Detector)

Catalog Number: 83682-3-PBS



Basic Information

Catalog Number:

83682-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10944

GenBank Accession Number:

BC009352

GeneID (NCBI):

55536

UNIPROT ID:

Q96GN5

Full Name:

cell division cycle associated 7-like

Calculated MW:

454 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

240639G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00669-1, CDCA7L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83682-1-PBS. Detection antibody: 83682-3-PBS. Standard: Ag10944. Range: 0.625-40 ng/mL.



Cytometric bead array standard curve of MP00669-2, CDCA7L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83682-2-PBS. Detection antibody: 83682-3-PBS. Standard: Ag10944. Range: 2.5-40 ng/mL.