

For Research Use Only

# Interferon alpha 1/IFNA1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83684-1-PBS



## Basic Information

Catalog Number:

83684-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_024013.3

GeneID (NCBI):

3447

UNIPROT ID:

P01562

Full Name:

interferon, alpha 13

Calculated MW:

22kDa

Purification Method:

Protein A purification

CloneNo.:

240667B6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00661-1, Interferon alpha 1/IFNA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83684-4-PBS. Detection antibody: 83684-1-PBS. Standard: Eg1177. Range: 0.156-20 ng/mL