

For Research Use Only

# CHI3L1/YKL40 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83685-1-PBS



## Basic Information

Catalog Number:

83685-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33123

GenBank Accession Number:

BC008568

GeneID (NCBI):

1116

UNIPROT ID:

P36222

Full Name:

chitinase 3-like 1 (cartilage  
glycoprotein-39)

Calculated MW:

383 aa, 43 kDa

Purification Method:

Protein A purification

CloneNo.:

240691C5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00663-1, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-1-PBS. Detection antibody: 83685-2-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00663-2, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-1-PBS. Detection antibody: 83685-3-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL