

For Research Use Only

CHI3L1/YKL40 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83685-3-PBS



Basic Information

Catalog Number:

83685-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33123

GenBank Accession Number:

BC008568

GeneID (NCBI):

1116

UNIPROT ID:

P36222

Full Name:

chitinase 3-like 1 (cartilage
glycoprotein-39)

Calculated MW:

383 aa, 43 kDa

Purification Method:

Protein A purification

CloneNo.:

240691E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00663-2, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-1-PBS. Detection antibody: 83685-3-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00663-3, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-2-PBS. Detection antibody: 83685-3-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL