

For Research Use Only

HSP47/SERPINH1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83689-2-PBS



Basic Information

Catalog Number:

83689-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001235.5

GeneID (NCBI):

871

UNIPROT ID:

P50454

Full Name:

serpin peptidase inhibitor, clade H
(heat shock protein 47), member 1,
(collagen binding protein 1)

Calculated MW:

46kDa

Purification Method:

Protein A purification

CloneNo.:

240752A12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00677-1, HSP47/SERPINH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83689-2-PBS. Detection antibody: 83689-1-PBS. Standard: Eg1201. Range: 0.391-50 ng/mL