

For Research Use Only

# RABGAP1L Recombinant antibody, PBS Only



Catalog Number: 83698-4-PBS

## Basic Information

|                                     |  |  |
|-------------------------------------|--|--|
| Catalog Number:<br>83698-4-PBS      | GenBank Accession Number:<br>BC041888              | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>9910                             | CloneNo.:<br>240663G12                         |
| Source:<br>Rabbit                   | UNIPROT ID:<br>B7ZAP0                              |  |
| Isotype:<br>IgG                     | Full Name:<br>RAB GTPase activating protein 1-like |  |
| Immunogen Catalog Number:<br>AG4856 | Calculated MW:<br>815 aa, 93 kDa                   |  |

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83698-4-RR (RABGAP1L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83698-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83698-4-RR against Human RABGAP1L were performed. The affinity constant is below 1 pM.