

For Research Use Only

GBP4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83699-3-PBS

Basic Information

Catalog Number:

83699-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12006

GenBank Accession Number:

BC070055

GeneID (NCBI):

115361

UNIPROT ID:

Q96PP9

Full Name:

guanylate binding protein 4

Calculated MW:

640 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240730A12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00655-1, GBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83699-1-PBS. Detection antibody: 83699-3-PBS. Standard: Ag12006. Range: 0.78-50 ng/mL.



Cytometric bead array standard curve of MP00655-2, GBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83699-2-PBS. Detection antibody: 83699-3-PBS. Standard: Ag12006. Range: 3.13-100 ng/mL.