

For Research Use Only

TNFSF15 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83700-1-PBS



Basic Information

Catalog Number:

83700-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG32092

GenBank Accession Number:

BC074941

GeneID (NCBI):

9966

UNIPROT ID:

O95150

Full Name:

tumor necrosis factor (ligand)
superfamily, member 15

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

240783A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00659-1, TNFSF15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83700-2-PBS. Detection antibody: 83700-1-PBS. Standard: Eg32092. Range: 0.781-100 ng/mL