

For Research Use Only

Mouse Haptoglobin Recombinant antibody, PBS Only (Capture)

Catalog Number:83702-2-PBS



Basic Information

Catalog Number: 83702-2-PBS	GenBank Accession Number: NM_017370.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 15439	CloneNo.: 240785C6
Source: Rabbit	UNIPROT ID: Q61646	
Isotype: IgG	Full Name: haptoglobin	
Immunogen Catalog Number: EG1197	Calculated MW: 39kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00665-2, MOUSE Haptoglobin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83702-2-PBS. Detection antibody: 83702-1-PBS. Standard: Eg1197. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00665-3, MOUSE Haptoglobin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83702-2-PBS. Detection antibody: 83702-4-PBS. Standard: Eg1197. Range: 1.25-40 ng/mL