

For Research Use Only

# GALNT2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83703-1-PBS

## Basic Information

Catalog Number:	83703-1-PBS	GenBank Accession Number:	BC041120	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	2590	CloneNo.:	240723D8
Source:	Rabbit	UNIPROT ID:	Q10471		
Isotype:	IgG	Full Name:	UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylglucosaminyltransferase 2 (GalNAc-T2)		
Immunogen Catalog Number:	AG11523	Calculated MW:	571 aa, 65 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00666-1, GALNT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83703-1-PBS. Detection antibody: 83703-2-PBS. Standard: Ag11523. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP00666-2, GALNT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83703-3-PBS. Detection antibody: 83703-1-PBS. Standard: Ag11523. Range: 6.25-100 ng/mL