

For Research Use Only

GALNT2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83703-3-PBS

Basic Information

Catalog Number:

83703-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11523

GenBank Accession Number:

BC041120

GeneID (NCBI):

2590

UNIPROT ID:

Q10471

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetyl(galactosaminyltransferase 2 (GalNAc-T2)

Calculated MW:

571 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

240723E6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00666-2, GALNT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83703-3-PBS. Detection antibody: 83703-1-PBS. Standard: Ag11523. Range: 6.25-100 ng/mL.