

For Research Use Only

ULBP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83709-2-PBS



Basic Information

Catalog Number:

83709-2-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0858

GenBank Accession Number:

NM_025217.4

GeneID (NCBI):

80328

UNIPROT ID:

Q9BZM5

Full Name:

UL16 binding protein 2

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

240413E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01920-1, ULBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83709-2-PBS. Detection antibody: 83709-4-PBS. Standard: Eg0858. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83709-2-PBS against Human ULBP2 was performed. The affinity constant is 2.30 nM.