

For Research Use Only

Acetylcholinesterase Recombinant antibody, PBS Only (Detector)

Catalog Number: 83714-2-PBS



Basic Information

Catalog Number:

83714-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12146

GenBank Accession Number:

BC094752

GeneID (NCBI):

43

UNIPROT ID:

P22303

Full Name:

acetylcholinesterase (Yt blood group)

Calculated MW:

614aa, 68 kDa; 640aa, 71 kDa

Purification Method:

Protein A purification

CloneNo.:

240640H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00685-1, ACHE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83714-1-PBS. Detection antibody: 83714-2-PBS. Standard: Ag12146. Range: 1.25-80 ng/mL.