

For Research Use Only

# MIB1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83718-2-PBS

## Basic Information

Catalog Number: 83718-2-PBS	GenBank Accession Number: BC022403	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57534	CloneNo.: 240676E6
Source: Rabbit	UNIPROT ID: Q86YT6	
Isotype: IgG	Full Name: mindbomb homolog 1 (Drosophila)	
Immunogen Catalog Number: AG2507	Calculated MW: 110 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00718-2, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-2-PBS. Detection antibody: 83718-4-PBS. Standard: Ag2507. Range: 1.25-80 ng/mL.