

For Research Use Only

TNF-alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 83721-2-PBS



Basic Information

Catalog Number:

83721-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000594.4

GeneID (NCBI):

7124

ENSEMBL Gene ID:

ENSG00000232810

UNIPROT ID:

P01375

Full Name:

tumor necrosis factor (TNF
superfamily, member 2)

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

240668H1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00702-1, TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83721-3-PBS. Detection antibody: 83721-2-PBS. Standard: Eg0816. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00702-3, TNF-alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83721-1-PBS. Detection antibody: 83721-2-PBS. Standard: Eg0816. Range: 0.628-40 ng/mL