

For Research Use Only

# IL-1RAP Recombinant antibody, PBS Only (Capture)

Catalog Number: 83726-4-PBS



## Basic Information

Catalog Number:

83726-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002182.4

GeneID (NCBI):

3556

UNIPROT ID:

Q9NPH3-1

Full Name:

interleukin 1 receptor accessory protein

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

240756H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00708-1, IL-1RAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83726-4-PBS. Detection antibody: 83726-2-PBS. Standard: Eg0827. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00708-2, IL-1RAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83726-4-PBS. Detection antibody: 83726-3-PBS. Standard: Eg0827. Range: 0.156-20 ng/mL