

For Research Use Only

IL1RAP Recombinant antibody, PBS Only (Capture)

Catalog Number: 83726-7-PBS



Basic Information

Catalog Number:

83726-7-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0827

GenBank Accession Number:

NM_002182.4

GeneID (NCBI):

3556

UNIPROT ID:

Q9NPH3-1

Full Name:

interleukin 1 receptor accessory protein

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

240756A12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00708-4, Human IL-1RAP Recombinant Matched Antibody Pair - PBS only. 83726-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0827. 83726-5-PBS was HRP conjugated as the detection antibody. Range: 7.8-500 pg/mL



Plasma of six individual healthy human donors was measured. The human IL-1RAP concentration of detected samples was determined to be 266.79 pg/mL with a range of 13.79 - 1,302.80 pg/mL.



Serum of six individual healthy human donors was measured. The human IL-1RAP concentration of detected samples was determined to be 21.67 pg/mL with a range of 7.89 - 42.37 pg/mL.