

For Research Use Only

# VEGF-C Recombinant antibody, PBS Only (Capture)

Catalog Number: 83727-2-PBS



## Basic Information

Catalog Number:

83727-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC035212

GeneID (NCBI):

7424

UNIPROT ID:

P49767

Full Name:

vascular endothelial growth factor C

Calculated MW:

419 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

240772F11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00710-1, VEGF-C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83727-2-PBS. Detection antibody: 83727-3-PBS. Standard: Eg32025. Range: 0.156-20 ng/mL