

For Research Use Only

# ACOX1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83731-3-PBS



## Basic Information

Catalog Number:

83731-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1405

GenBank Accession Number:

BC010425

GeneID (NCBI):

51

UNIPROT ID:

Q15067

Full Name:

acyl-Coenzyme A oxidase 1,  
palmitoyl

Calculated MW:

74 kDa

Purification Method:

Protein A purification

CloneNo.:

240754D4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00709-1, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-4-PBS. Detection antibody: 83731-3-PBS. Standard: Ag1405. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00709-2, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-2-PBS. Detection antibody: 83731-3-PBS. Standard: Ag1405. Range: 1.56-100 ng/mL