

For Research Use Only

ARID1A Recombinant antibody, PBS Only (Detector)

Catalog Number: 83733-4-PBS



Basic Information

Catalog Number:

83733-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30749

GenBank Accession Number:

NM_006015

GeneID (NCBI):

8289

UNIPROT ID:

O14497

Full Name:

AT rich interactive domain 1A (SWI-like)

Calculated MW:

242 kDa

Purification Method:

Protein A purification

CloneNo.:

240690G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00714-1, ARID1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83733-1-PBS. Detection antibody: 83733-4-PBS. Standard: Ag30749. Range: 0.625 -80 ng/mL