

For Research Use Only

# USP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83734-1-PBS



## Basic Information

Catalog Number:

83734-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28297

GenBank Accession Number:

BC125131

GeneID (NCBI):

7375

UNIPROT ID:

Q13107

Full Name:

ubiquitin specific peptidase 4 (proto-oncogene)

Calculated MW:

963 aa, 109 kDa

Purification Method:

Protein A purification

CloneNo.:

240797A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00717-2, USP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83734-3-PBS. Detection antibody: 83734-1-PBS. Standard: Ag28297. Range: 0.625-80 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83734-1-PBS against Human USP4 was performed. The affinity constant is 2.69 nM.