

For Research Use Only

# Osteoprotegerin/TNFRSF11B Recombinant antibody, PBS Only (Capture)

Catalog Number: 83737-1-PBS



## Basic Information

Catalog Number:

83737-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_002546.4

GeneID (NCBI):

4982

UNIPROT ID:

O00300

Full Name:

tumor necrosis factor receptor  
superfamily, member 11b

Calculated MW:

46kDa

Purification Method:

Protein A purification

CloneNo.:

240671A4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00725-2, Osteoprotegerin/TNFRSF 11B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83737-1-PBS. Detection antibody: 83737-2-PBS. Standard: Eg0976. Range: 0.156-20 ng/mL