

For Research Use Only

Osteoprotegerin/TNFRSF11B Recombinant antibody, PBS Only (Detector)

Catalog Number: 83737-4-PBS



Basic Information

Catalog Number:

83737-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002546.4

GeneID (NCBI):

4982

UNIPROT ID:

O00300

Full Name:

tumor necrosis factor receptor
superfamily, member 11b

Calculated MW:

46kDa

Purification Method:

Protein A purification

CloneNo.:

240671A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00725-1, Osteoprotegerin/TNFRSF 11B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83737-3-PBS. Detection antibody: 83737-4-PBS. Standard: Eg0976. Range: 0.156-20 ng/mL