

For Research Use Only

# Mouse CX3CL1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83738-1-PBS

## Basic Information

Catalog Number: 83738-1-PBS	GenBank Accession Number: NM_009142.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 20312	CloneNo.: 240738A3
Source: Rabbit	UNIPROT ID: O35188	
Isotype: IgG	Full Name: chemokine (C-X3-C motif) ligand 1	
	Calculated MW: 42 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
mouse

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00688-2, MOUSE CX3CL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83738-1-PBS. Detection antibody: 83738-3-PBS. Standard: Eg0704. Range: 1.25-80 ng/mL